


# Utilization of Tumor-derived Extracellular Vesicles for Patient Stratification and Biomarker Discovery

Claire Seibold, Tia Seibold, Margie Kinnersley, Tre Blohm, Kelley Van Vaerenberghe, Rachel Short-Miller, Sean Lodmell, Amanda Mast, Riley Kemp, and Katie Havranek



## Introduction



- Liquid biopsy is poised to play an increasingly important role in clinical decision making by enabling earlier cancer detection, diagnosis, patient stratification, and treatment monitoring.
- Most approaches rely on detection of cell-free circulating tumor DNA or circulating tumor cells, but limitations exist, including minuscule abundance, unfavorable signal to noise ratios, and limited insight into mechanisms driving disease.
- Extracellular vesicles (EVs) are an alternative liquid biopsy analyte comprising lipid membrane encapsulated particles carrying diverse protein, nucleic acid, and metabolite cargos.
- We have developed a clinically applicable, multi-omic EV subpopulation interrogation pipeline that robustly profiles tumor and brain derived EVs (TDEVs and BDEVs) in biofluids utilizing FYR's novel EV enrichment technology called SPARCs.
- We have applied our EV-Omic (EVO) pipeline in combination with machine learning to distinguish adult high-grade gliomas and benign brain tumors from healthy patients, using blood-based EV proteomic and transcriptomic signatures.



## Methods

Plasma was processed from whole blood in Streck Cell-Free BCT tubes (Streck, LaVista, NE). Total EVs were isolated via ion-exchange chromatography and characterized in concordance with the Minimal Information for Studies of Extracellular Vesicles (MISEV) 2023 guidelines. EVs were incubated with Tumor SPARCs and Neuro SPARCs to enrich tumor-derived EVs and brain-derived EVs, respectively. EV subpopulations were split for RNA and protein downstream processing. Small RNA was sequenced on an Element Biosciences AVITI (San Diego, CA). Protein was subjected to LC-MS/MS on an Orbitrap Astral (Thermo Scientific Fisher, Waltham, MA) at Cedars Sinai Precision Biomarker Laboratories (Beverly Hills, CA). Machine Learning model training utilized 10-fold cross validation, and the best model was chosen based on AUC. Candidate biomarker features were chosen from the leading model AUC.

## Study Design and EV Characterization

	N	Specimen Providers	Age				Sex		Race/Ethnicity		
			<51	51-60	61-70	71-80	Male	Female	White	Asian	Black
<b>Glioma</b>	<b>52</b>	2	18	14	18	2	23	29	51	1	0
<b>Benign</b>	<b>15</b>	1	6	3	5	1	2	13	15	0	0
<b>Healthy</b>	<b>38</b>	2	13	12	11	2	21	17	34	0	4



Table 1. Patients included in the study.



EV marker protein composition was confirmed via Capillary Electrophoresis (CE) Western Blot. EV markers CD63 (MISEV Cat1a),  $\beta$ -Actin (MISEV Cat2b), CD41a (MISEV Cat1b), and Albumin (Cat3a), confirm the presence of EVs as well as the common plasma co-isolate Albumin.

Nanoparticle Tracking Analysis (NTA) yielded an average of  $5.7 \times 10^{11}$  particles/mL plasma with a mean diameter of 145.7nm. There was a significant difference in concentration between Brain Cancer and Healthy donors where  $p=0.036$  via Kruskal-Wallis test. This is consistent with literature citing elevated circulating EV content in brain cancer patient plasma.

## Brain Tumor Plasma Proteomic Signatures



Differential expression analysis of brain tumor (malignant and benign) patients relative to healthy controls identified 256 significantly enriched proteins and 154 significantly depleted proteins using Tumor SPARCs and 598 significantly enriched proteins and 360 significantly depleted proteins using Neuro SPARCs (Fig.3). Comparison of identified DEPs from Tumor SPARCs vs Neuro SPARCs demonstrates that both panels provide unique DEP information (Fig.4). Neuro SPARCs enrichment demonstrated increased identification of CNS proteins (classified as brain elevated according to the Human Protein Atlas) (Fig.5). Heatmap hierarchical clustering by disease status of Tumor SPARCs DEPs and Neuro SPARCs DEPs (Fig.6).

## Benign Tumor Discrimination






Figure 8. Gene ontology terms listing the top 5 most significant terms for each GO category (cellular function, biological process) for DEPs identified in differential expression analysis.

## Biomarker Discovery


Objective: Distinction of brain tumors (malignant and benign) from healthy patients.

	Protein Tumor SPARCs	Protein Neuro SPARCs	miRNA Tumor SPARCs	miRNA Neuro SPARCs
<b>AUC</b>	0.98	0.96	0.80	0.83
<b>Accuracy</b>	0.87	0.88	0.77	0.74
<b>Sensitivity</b>	0.82	0.82	0.47	0.56
<b>Specificity</b>	0.95	0.97	0.89	0.89
<b>Features</b>	16	17	13	17




Objective: Distinction of glioblastomas from healthy patients.

	Protein Tumor SPARCs	Protein Neuro SPARCs	miRNA Tumor SPARCs	miRNA Neuro SPARCs
<b>AUC</b>	0.93	0.94	0.86	0.83
<b>Accuracy</b>	0.89	0.84	0.79	0.74
<b>Sensitivity</b>	0.75	0.72	0.78	0.56
<b>Specificity</b>	1	0.95	0.79	0.89
<b>Features</b>	11	16	17	18




## Subtyping Pilot: Tissue MGMT Status






For n=6 glioblastoma patient plasma specimens, matched tissue was assayed for MGMT promoter methylation status using PCR. Neuro SPARCs proteomic profiles provide preliminary evidence supporting the ability of BDEV subpopulations to identify novel liquid biomarkers for MGMT methylation status.

## Longitudinal Pilot: Glioblastoma



	Age	Timepoints 1, 2, 3	Treatment	Response
<b>Donor 127</b>	46	✓, -, -	NA	NA
<b>Donor 134</b>	58	✓, -, ✓	Radiation	Stable
<b>Donor 140</b>	66	✓, ✓, ✓	Radiation + Temozolomide	Stable
<b>Donor 198</b>	46	✓, ✓, ✓	Radiation + Temozolomide	Progressing (leukopenia, anemia)
<b>Donor 284</b>	69	✓, ✓, -	NA	NA



## Conclusions and Future Directions

- EV subpopulations in brain tumor patient plasma provide opportunities to identify disease via minimally invasive methods and provide a better understanding of tumor biology.
- Tumor SPARCs and Neuro SPARCs capture EV subpopulations with unique characteristics, but Tumor SPARCs information alone provided slightly improved performance scores. Combining scores and markers from both panel models is under investigation for ability to improve patient stratification and biomarker discovery.
- Preliminary data indicate that SPARCs will prove useful in blood-based tissue subtyping, as well as longitudinal monitoring of patients.

